

Laboratory Name:	
Laboratory Address:	
Date of this Packet:	
Insert Revision:	05/27/2014

Fisher Scientific Sure-Vue® Strep A Test (27 tests) No. 23900532 Laboratory Procedure

This procedure is intended to provide a ready outline reference for performance of the assay. These abbreviated directions for use are not intended to replace the complete package insert. **Any modifications to this document are the sole responsibility of the Facility.**

This is a Waived Complexity test.

1. Intended Use

The **Sure-Vue® Strep A Test** is a rapid chromatographic immunoassay for the qualitative detection of Strep A antigen from throat swab specimens to aid in the diagnosis of Group A Streptococcal infection.

2. Summary

Streptococcus pyogenes is non-motile gram-positive cocci, which contains the Lancefield group A antigen that can cause serious infections such as pharyngitis, respiratory infection, impetigo, endocarditis, meningitis, puerperal sepsis, and arthritis.¹ Left untreated, these infections can lead to serious complications, including rheumatic fever and peritonsillar abscess.² Traditional identification procedures for Group A Streptococci infection involve the isolation and identification of viable organisms using techniques that require 24 to 48 hours or longer.³

The **Sure-Vue® Strep A Dipstick** is a rapid test to qualitatively detect the presence of Strep A antigen in throat swab specimens, providing results within 5 minutes. The test utilizes antibodies specific for whole cell Lancefield Group A *Streptococcus* to selectively detect Strep A antigen in a throat swab specimen.

3. Test Principle

The **Sure-Vue® Strep A Test** is a qualitative, lateral flow immunoassay for the detection of Strep A carbohydrate antigen in a throat swab. In this test, antibody specific to Strep A carbohydrate antigen is coated on the test line region of the strip. During testing, the extracted throat swab specimen reacts with an antibody to Strep A that is coated onto particles. The mixture migrates up the membrane to react with the antibody to Strep A on the membrane and generates a red line in the test region. The presence of this red line in the test region indicates a positive result, while its absence indicates a negative result. To serve as a procedural control, a red line will always appear in the control region if the test has been performed properly. If a red control line does not appear, the test result is not valid.

4. Specimen Collection and Handling

Specimen:	Collect the throat swab specimen with the sterile swab that is provided in the kit. Rayon transport swabs containing modified Stuart's or Amies liquid medium can also be used with this product. Swab posterior pharynx, tonsils and other inflamed areas. Avoid touching the tongue, cheeks and teeth with the swab ⁴ .
Specimen Storage:	Testing should ideally be performed immediately after the specimens have been collected. Swab specimens may be stored in a clean, dry plastic tube for up

	to 8 hours at room temperature or 72 hours at 2-8°C (36-46°F).
Handling Precautions:	Only use reagents provided in the kit. If a culture is desired, lightly roll the swab tip onto a Group A selective blood agar plate before using the swab in the Sure-Vue® Strep A Test . CAUTION: Handle all specimens as if they contain infectious agents.

5. Reagents and Materials

A. Materials Provided

Component	Content	Quantity
Test strips		27
Sterile swabs		27
Disposable extraction test tubes		27
Reagent A	2M Sodium Nitrite	10 mL
Reagent B	0.2M Acetic Acid	10 mL
Positive control	Non-viable Strep A, 0.09% NaN ₃	1 mL
Negative control	Non-viable Strep C, 0.09% NaN ₃	1 mL
Workstation		1
Package insert		1

B. Materials Required But Not Provided

- Timer

C. Storage and Stability

The kit can be stored at room temperature or refrigerated (2-30°C). The test strip must remain in the sealed pouch until use. **DO NOT FREEZE**. The test strip and the reagents are stable through the expiration date printed on the box. Do not use beyond the expiration date.

6. Quality Control

Internal Quality Control

Internal procedural controls are included in the test. A red line appearing in the control region (C) is an internal positive procedural control. It confirms sufficient specimen volume and correct procedural technique. A clear background is an internal negative background control. If the test is working properly, the background in the result area should be white to light pink and not interfere with the ability to read the test result.

External Quality Control

It is recommended that a positive and negative external control be run once per kit, and as deemed necessary by your internal laboratory procedures. External positive and negative controls are supplied in the kit. Alternatively, other Group A and non-Group A *Streptococcus* ATCC reference strains may be used as external controls. Some commercial controls may contain interfering preservatives; therefore, other commercial controls are not recommended.

Procedure for External Quality Control Testing

1. Add 4 full drops of Reagent A and 4 full drops of Reagent B into an extraction test tube. Tap the bottom of the tube gently to mix the liquid.
2. Add 1 full drop of positive or negative control solution into the tube, holding the bottle upright.
3. Place a clean swab into the tube. Rotate the swab 10 times in the tube. Leave the swab in the tube for 1 minute. Then press the swab against the side of the tube and squeeze the bottom of the tube while removing the swab so that most of the liquid stays in the tube. Discard the swab.
4. Continue with Step 4 of the Test Procedure.

7. Precautions


- For professional *in vitro* diagnostic use only. Do not use after expiration date.
- Do not eat, drink or smoke in the area where the specimens and kits are handled.

- Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout the procedure and follow the standard procedures for proper disposal of specimens.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
- **WARNING:** Reagent A is harmful if swallowed or absorbed through skin. May cause eye irritation.
- **CAUTION:** Reagent B may cause skin, eye and respiratory tract irritation.
- The positive and negative controls contain sodium azide (NaN₃) as a preservative.
- Do not interchange reagent bottle caps.
- Do not interchange external control solution bottle caps.
- Humidity and temperature can adversely affect results.

8. Test Procedure

Allow the test strip, reagents, and/or controls to reach room temperature (15-30°C) prior to testing.

1. Remove the test strip from the sealed foil pouch and use it as soon as possible. Best results will be obtained if the test is performed immediately after opening the foil pouch.
2. Hold the Reagent A bottle upright and add 4 full drops (approximately 240 µL) to an extraction test tube. Reagent A is red in color. Hold the Reagent B bottle upright and add 4 full drops (approximately 160 µL) to the tube. Reagent B is colorless. The addition of Reagent B to Reagent A changes the color of the solution from red to pale yellow. Tap the bottom of the tube gently to mix the liquid.
3. Immediately add the throat swab into the tube of pale yellow solution. Rotate the swab vigorously 10 times in the tube. Leave the swab in the tube for 1 minute. Then press the swab against the side of the tube and squeeze the bottom of the tube while removing the swab so that most of the liquid stays in the tube. Discard the swab.
4. With arrows pointing down, place the test strip into the tube of solution and then start the timer. If the procedure is followed down, correctly, the liquid should be at or just below the maximum line (MAX) on the test strip. See the illustration below.
5. Leave the strip in the tube and read the result at 5 minutes. The result is invalid after 10 minutes.



8. Result Interpretation

(Please refer to the illustration above.)

POSITIVE*: **Two distinct red lines appear.** One line should be in the control region (C) and another line should be in the test region (T). A positive result indicates that Strep A was detected in the sample.

***NOTE:** The intensity of the red color in the test line region (T) will vary depending on the concentration of Strep A present in the sample. Therefore, any shade of red in the test region (T) should be considered positive.

NEGATIVE: One red line appears in the control region (C). No apparent red or pink line appears in the test region (T). A negative result indicates that Strep A is not present in the sample, or is present below the detectable level of the test. The patient's sample should be cultured to confirm the absence of Strep A infection. If clinical symptoms are not consistent with results, obtain another sample for culture.

INVALID: Control line fails to appear. Insufficient sample volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test strip. If the problem persists, discontinue using the test kit immediately and call 1-866-216-0094 for Technical Assistance.

9. Limitations

1. The **Sure-Vue® Strep A Test** is for *in vitro* diagnostic use only. The test should be used for the detection of Strep A antigen in throat swab specimens only. Neither the quantitative value nor the rate of increase in Strep A antigen concentration can be determined by this qualitative test.
2. This test will only indicate the presence of Strep A antigen in the specimen from both viable and non-viable Group A *Streptococcus* bacteria.
3. A negative result obtained from this kit should be confirmed by culture. A negative result may be obtained if the concentration of the Strep A antigen present in the throat swab is not adequate or is below the detectable level of the test.
4. Excess blood or mucus on the swab specimen may interfere with test performance and may yield a false positive result. Avoid touching the tongue, cheeks, and teeth⁴ and any bleeding areas of the mouth with the swab when collecting specimens.
5. As with all diagnostic tests, all results must be interpreted together with other clinical information available to the physician.

10. Expected Values

Approximately 15% of pharyngitis in children ages 3 months to 5 years is caused by Group A beta-hemolytic *Streptococcus*.⁵ In school-aged children and adults, the incidence of Strep throat infection is about 40%.⁶ This disease usually occurs in the winter and early spring in temperate climates.³

11. Performance Characteristics

Using three medical centers for evaluation, a total of 499 throat swabs were collected from patients exhibiting symptoms of pharyngitis. Each swab was rolled onto a sheep blood agar plate, and then tested by the **Sure-Vue® Strep A Test**. The plates were further streaked for isolation, and then incubated at 37°C with 5-10% CO₂ and a Bacitracin disk for 18-24 hours. The negative culture plates were incubated for an additional 18-24 hours. Possible GAS colonies were subcultured and confirmed with a commercially available latex agglutination grouping kit.

Of the 499 total specimens, 375 were found to be negative by culture and 124 were found to be positive by culture. During this study, two Strep F specimens yielded positive results with the Test. One of these specimens was re-cultured, then re-tested and yielded a negative result. Three additional different Strep F strains were cultured and tested for cross-reactivity and also yielded negative results.

		Culture	
		+	-
Sure-Vue® Strep A	+	120	20

Test	-	4	355
-------------	---	---	-----

Sensitivity: 120/124 = 97% (91% to 99%)*
Specificity: 355/375 = 95% (92% to 97%)*
Accuracy: 475/499 = 95% (93% to 97%)*
Prevalence: 124/499 = 25%
PPV (+): 120/140 = 86% (79% to 91%)*
NPV (-): 355/359 = 99% (97% to 100%)*
 * Denotes a 95% Confidence Interval

Positive Culture Classification	Sure-Vue® Strep A Test/Culture	% Correct
Rare	10/11	91%
1+	9/9	100%
2+	17/19	89%
3+	36/37	97%
4+	48/48	100%

12. Cross-Reactivity

The following organisms were tested at 1.0×10^7 organisms per test and were all found to be negative when tested with the **Sure-Vue® Strep A Test**. No mucoid-producing strains were tested.

Group B <i>Streptococcus</i>	Group C <i>Streptococcus</i>
Group F <i>Streptococcus</i>	Group G <i>Streptococcus</i>
<i>Streptococcus pneumoniae</i>	<i>Streptococcus sanguis</i>
<i>Streptococcus mutans</i>	<i>Enterococcus faecalis</i>
<i>Staphylococcus aureus</i>	<i>Staphylococcus epidermidis</i>
<i>Corynebacterium diphtheriae</i>	<i>Serratia marcescens</i>
<i>Candida albicans</i>	<i>Klebsiella pneumoniae</i>
<i>Pseudomonas aeruginosa</i>	<i>Bordetella pertussis</i>
<i>Neisseria meningitidis</i>	<i>Neisseria gonorrhoeae</i>
<i>Neisseria sicca</i>	<i>Neisseria subflava</i>
<i>Branhamella catarrhalis</i>	<i>Haemophilus influenza</i>

13. POL Studies

Three physicians' offices were used to conduct an evaluation of the **Sure-Vue® Strep A Test**. Personnel with various educational backgrounds performed the testing. Each physician's office tested a randomly coded panel of samples consisting of negative (20), low positive (20), and medium positive (20) for three days. The results obtained had a 96% correlation with the expected results.

14. References

- Manual of Clinical Microbiology, 6th Edition, ASM Press, p. 299-307.
- Webb, KH. *Pediatrics* (Feb 1998), 101: 2, 2.
- Bisno AL, Gerber MA, Gwaltney JM, Kaplan EL, Schwartz RH. *Clinical Infectious Diseases* (1997), 25, 574-83.
- Shea, Y.R., Specimen Collection and Transport, in *Clinical Microbiology Procedures Handbook*, Isenberg, H.D., American Society of Microbiology, Washington, D.C., 1.1.1-1.1.30, 1992.
- Nussinovitch, M, Finkelstein Y, Amir J, Varsano, I. *Clinical Pediatrics* (June 1999), 357-360.
- Woods WA, Carter CT, Stack M, Connors Jr AF, Schlager TA, *Southern Medical Journal* (May 1999), 491-492.

▪

Test Procedure Approval and Review Sheet

Prepared By:	
Date:	
Supervisor Review:	
Date:	
Laboratory Director or Designee Approval:	
Implementation Date:	
Supersedes Procedure Dated:	
Date Procedure Retired:	

Laboratory Director or Designee	Date Reviewed	Laboratory Director or Designee	Date Reviewed

Sure-Vue® Strep A Test Method Verification Form

Account Name: _____

Address: _____

Telephone: _____

**Sure-Vue™ Strep A
Test Lot #/Exp:** _____

Date: _____

Supervisor Signature: _____

Record the results from reference samples below.

Record the Sample #, the **Sure-Vue® Strep A Test** results, Tester's Initials, and any comments. After the **Sure-Vue® Strep A Test** results have been recorded (positive or negative) then record the Expected Results (positive or negative).

Sample #	Expected Results	Sure-Vue® Strep A Test Result	Tester's Initials	Comments

Sure-Vue® Strep A Test Method Verification Form (continued)

Sample #	Expected Results	Sure-Vue® Strep A Test Result	Tester's Initials	Comments

Review: _____ Date: _____

Laboratory Director Review and Approval for Clinical Use: _____

Date: _____

Sure-Vue® Strep A Test External Quality Control

Name of Facility: _____

External QC testing is recommended:

- Be run once per kit, and as deemed necessary by your internal laboratory procedures.
- When required by local, state, and/or federal regulations, accrediting groups, or your lab's Quality Control procedures.

Date	Sure-Vue® Strep A Test Kit Lot/Exp	Positive Ctrl Lot/Exp	Negative Ctrl Lot/Exp	Positive Result	Negative Result	Tester's Initials	Comments

Reviewed by: _____

Date: _____

Sure-Vue® Strep A Test Quality Control and Patient Record

Lot Number _____ Exp. Date _____

Recommendation is that external positive and negative controls be run with each new kit.

Record the Date, Patient's Name, Patient Test Result, Internal Control Results and the performer's initials.

Positive Internal Control=the red line appearing at the "control line" position; Negative Internal Control=background color should be white to light pink.

Date	Patient Name	Patient ID Number	Patient Results	Are Internal Control Results Invalid or Valid?		Internal Control Results		Comments	Tech
				Invalid	Valid	+	-		

Reviewed By: _____

Date: _____

Quality Assessment Review Form and Checklist

These forms are used for periodical review of the patient testing process. These should be filed with the quality assessment records.

Quality Assessment Activity Initial	Comments	Date	Initials
Patient Test Management: Evaluate criteria for specimen submission, handling, and rejection; test results requisitions and reporting, accuracy and reliability of reports.			
Quality Control: Assess calibration and control data, reference range verification, errors in reporting results, corrective actions taken with appropriate documentation records.			
Proficiency Testing: Review the effectiveness of corrective actions taken for unsatisfactory performance or failures.			
Comparison of Test Results: Review at least semi-annually comparative results for multiple methods, instruments, or site correlations when more than one procedure exists.			
Relationship of Patient Test Information to Test Results: Evaluate patient test reports for accuracy of patient information, test results, and normal ranges. Identify and evaluate results inconsistent with Patient's age, sex, diagnosis, and other test parameters.			
Personnel: Evaluate the effectiveness of policies and procedures for assuring employees competence of testing and reporting test results.			
Communications: Evaluate documented problems and corrective actions that occur between the laboratory and the authorized individual who orders or receives the test result.			
Complaint Investigation: Evaluate documented complaints and corrective actions.			
Quality Assessment Reviews with Staff: Document discussion with Staff regarding identified problems and corrective actions during the QA review.			

Corrective Action Form

Problem/Error

Corrective Action

Problem/Error	Corrective Action

Technologist: _____

Date: _____

Supervisor: _____

Date: _____

Laboratory Director: _____

Date: _____

TEMPERATURE LOG

Equipment: _____

Name of Facility: _____

To be recorded at the beginning of each workday. Temperature Range: _____

Date	°C	Initials	Adjustments	Date	°C	Initials	Adjustments

Tips for Successful PT Performance

- Strictly follow the PT provider's storage or handling requirement ***before testing PT specimens***.
- Analyze PT specimens ***within the time frame*** provided by the PT provider.
- Contact the PT provider ***promptly*** when specimens are received damaged. You may be able to receive a replacement immediately.
- Avoid clerical error when filling out PT answer sheets. Be sure to ***enter the correct result next to the correct analyte*** on the answer form.
- Remember to identify the instrument or method you are using to perform your PT so you are ***graded among your peer group***.
- Make copies of all answer forms ***before submitting them*** to your PT provider.
- Please contact Technical Support at 877-441-7440 or Lateral.Flow.Support@alere.com for further information on proficiency providers.

Certification of Training

This is to verify that personnel responsible for running the **Sure-Vue® Strep A Test** at _____ have been thoroughly in-serviced on the test and the test procedure. This has included:

- **Review of the package insert**
- **Demonstration of the product assay**
- **Successful performance of the Sure-Vue® Strep A Test and interpretation of results**

Names of the personnel who have been trained with the **Sure-Vue® Strep A Test** and are responsible for reporting patient results:

PRINT NAME	SIGNATURE	DATE

Signature of Laboratory Director(s) responsible for personnel and testing:

Signature

Date

Signature

Date

Trainer

Date

Testing Personnel Training Assessment

Test Method: Sure-Vue® Strep A Test

Procedure	Satisfactory	Unsatisfactory	Not Applicable	Comments / Corrective Actions
<i>Observation of Test Performance:</i>				
Patient Sample Preparation (if applicable)				
Specimen Handling/Processing				
Testing				
Recording/Reporting Results				
<i>Assessment of Test Performance Using Known Samples</i>				
<i>Review of Records:</i>				
Patient/Quality Control Log Sheet Records				
Proficiency Testing Records				
<i>Assessment of Problem Solving Skills</i>				

(Attach all supporting documents)

Evaluator: _____

Date: _____

Employee: _____

Sure-Vue® Strep A Test Quiz

Name: _____

Date: _____

Circle T (True) or F (False) for each Question:

- | | | |
|--|---|---|
| 1. The Sure-Vue® Strep A Test must be refrigerated at 2-8°C. | T | F |
| 2. The Sure-Vue® Strep A Test strip foil pouch may be left open for extended periods of time without affecting the test results. | T | F |
| 3. Throat swab is the only acceptable sample type for testing on the Sure-Vue® Strep A Test . | T | F |
| 4. Throat swab specimens may be held in a clean dry plastic tube up to 72 hours at 2-8°C before performing the Sure-Vue® Strep A Test . | T | F |
| 5. The Sure-Vue® Strep A Test does not require live organism to be positive. | T | F |
| 6. It is recommended that the swab provided in the kit be used. | T | F |
| 7. Testing should ideally be performed 10 minutes after the specimens have been collected. | T | F |
| 8. Test Results should be read at 5 minutes. | T | F |
| 9. The appearance of a red Control line and a red Test line is a positive result. | T | F |
| 10. External Positive and Negative Controls should be run with each new lot, each new shipment, monthly as a check on storage, each untrained operator and as otherwise required by your lab internal quality system procedures. | T | F |

Sure-Vue® Strep A Test Quiz Answers

	Answer Key	Explanation
1.	F	The Sure-Vue® Strep A Test can be stored at room temperature or refrigerated: 2°-30°C (36°-86°F).
2.	F	Immediately close the Sure-Vue® Strep A Test strip foil pouch after removing a test strip.
3.	T	Throat swab specimen is the only acceptable sample type for the Sure-Vue® Strep A Test .
4.	T	Throat swab specimens may be stored at room temperature for up to 8 hours or refrigerated at 2°-8°C for up to 72 hours prior to testing.
5.	T	The Sure-Vue® Strep A Test does not require live organism to be positive.
6.	T	Collect the throat swab specimen with the sterile swab that is provided in the kit.
7.	F	Testing should ideally be performed immediately after the specimens have been collected.
8.	T	Read test results at 5 minutes. Do not read after 10 minutes.
9.	T	Two distinct red lines appear; one in the control region and another in the test region. Any shade of red in the test region should be considered positive.
10.	T	External Positive and Negative Controls should be run with each new lot, each new shipment, monthly as a check on storage, each untrained operator and as otherwise required by your lab internal quality system procedures.