

**USP
797**

DON'T TAKE REGULATORY RISKS

Make sure you comply with new manufacturing standards.

CONTAMINATION COSTS

1,084

adverse events
and **87**
deaths
from 2001-13¹

CONTAMINATION OF STERILE PRODUCTS

Most common compounding error²



People
#1



SOURCE


of cleanroom
contamination³



30% of microbial
detection in sterile
cleanrooms is on clothing⁴



>\$3.1MM
the cost of
contamination¹




The annual cost of
contamination at a single site⁵

PROTECT YOUR PROCESS



Open vials of injectable
pharmaceutical product are
at risk from contamination



FDA recommends improved
arm coverage to reduce
contamination risk⁶

STERILE, SINGLE-USE GARMENTS

deliver enhanced
protection, compliance
and optimal, predictable
performance

THE ASEPTIC SOLUTION

Kimtech Pure* A5 Sterile Cleanroom Apparel.

Cleanroom-level protection for compounding pharmacies.



✓ Coveralls

✓ Sleeves

✓ Hoods

✓ Masks

✓ Boots

✓ Gloves

✓ Shoe Covers

✓ And more...

Cleanroom Expertise for Compounding Pharmacies. Kimtech Pure* A5 Sterile Apparel

For more information, visit: www.kimtech.com

1. U.S. Illnesses and Deaths Associated with Compounded Medications, The Pew Charitable Trusts Health Initiatives, September 6, 2013. <http://www.pewhealth.org/other-resource/us-illnesses-and-deaths-associated-with-compounded-medications-85899468587>
2. Ibid
3. PDA Midwest Meeting 032510
4. PDA Midwest Meeting 2012, www.pda.org
5. Ibid
6. www.fda.org

THE ASEPTIC SOLUTION



Kimtech Pure® A5 Sterile Cleanroom Apparel**

Fisher Scientific Code	KCP Code	Description	Total/Case	Product Features
19-057-475	88800	Small Coverall, Sterile	25	<ul style="list-style-type: none"> • Cloth like comfort • Built in thumb loop • Patent-Pending Snap Technology
19-057-476	88801	Medium Coverall, Sterile	25	
19-057-490	88802	Large Coverall, Sterile	25	
19-057-491	88803	XL Coverall, Sterile	25	
19-057-492	88804	2XL Coverall, Sterile	25	
19-057-493	88805	3XL Coverall, Sterile	25	
19-057-494	88806	4XL Coverall, Sterile	25	
18-042-911	12917	5XL Coverall, Sterile	25	
18-042-910	12914	6XL – 8XL Coverall, Sterile	25	



19-048-606	36072	Universal Integrated Hood and Mask, Sterile	75	<ul style="list-style-type: none"> • Stretch fit elastic • Easier donning – 2 steps become 1 • Pull-down ties assist gowning process
------------	-------	---	----	---



17-150-200	36077	Sleeve Protector, Sterile	100 Pairs	<ul style="list-style-type: none"> • Blue-line donning signal feature • Thumb loop • 18" length
------------	-------	---------------------------	-----------	--



19-169-903	31696	Universal Boots w/ Grasp Ties and Vinyl Edge, Sterile	100 Pairs	<ul style="list-style-type: none"> • Side grasp ties for secure fit • Vinyl edge reduces slipping/liquid penetration
19-169-902	31683	SM/MD Boots w/ Grasp Ties and Vinyl Edge, Sterile	100 Pairs	
19-169-904	31697	XL/2XL Boots w/ Grasp Ties and Vinyl Edge, Sterile	100 Pairs	

** Gamma Irradiated SAL:10⁻⁶

Kimtech Pure® A8 Unitrax® Shoe Covers



Fisher Scientific Code	KCP Code	Description	Size	Color	Total/Case	Product Features
19-169-650	39370	A8 Shoe Covers w/ seamless bottom	Small/Med	White	300	<ul style="list-style-type: none"> • Dry particle and fluid protection • Extra grip • Low linting, Category 1 Helmke Drum Test
19-169-651	39371	A8 Shoe Covers w/ seamless bottom	Univ	White	300	
19-169-652	39372	A8 Shoe Covers w/ seamless bottom	XL/2XL	White	300	



KleenGuard® A20 Shoe Covers

Fisher Scientific Code	KCP Code	Description	Size	Color	Total/Case	Product Features
19-005-059	66857	A20 Shoe Covers w/ seamless elastic sole, toe & heel	XL/2XL	Blue	300	<ul style="list-style-type: none"> • Dry particle protection • No tack • Packed for cleanroom • Recommended use – under boot covers



KleenGuard® Bouffant Caps

Fisher Scientific Code	KCP Code	Description	Size	Color	Total/Case	Product Features
17-981-43A	36850	A10 Bouffant Cap	Medium - 21"	White	1,000	<ul style="list-style-type: none"> • Dry particle protection • Latex free • Packed for cleanroom
17-981-43B	36860	A10 Bouffant Cap	Large - 24"	White	1,000	
N/A	66829	A20 Bouffant Cap	Large - 24"	White	500	



Kimtech Pure® G3 Sterile Sterling® Nitrile Gloves

Packaged in a Class 100 Cleanroom

Fisher Scientific Code	KCP Code	Description	Size	Total/Case	Product Features
19-047-730	11821	Kimtech Pure® G3 Sterile Sterling® Nitrile Gloves - 12"	6	300	<ul style="list-style-type: none"> • Latex free • Static dissipative in use • Beaded cuff • Hand specific pairs
19-047-732	11822	Kimtech Pure® G3 Sterile Sterling® Nitrile Gloves - 12"	6.5	300	
19-047-733	11823	Kimtech Pure® G3 Sterile Sterling® Nitrile Gloves - 12"	7	300	
19-047-734	11824	Kimtech Pure® G3 Sterile Sterling® Nitrile Gloves - 12"	7.5	300	
19-047-735	11825	Kimtech Pure® G3 Sterile Sterling® Nitrile Gloves - 12"	8	300	
19-047-736	11826	Kimtech Pure® G3 Sterile Sterling® Nitrile Gloves - 12"	8.5	300	
19-047-737	11827	Kimtech Pure® G3 Sterile Sterling® Nitrile Gloves - 12"	9	300	
19-047-738	11828	Kimtech Pure® G3 Sterile Sterling® Nitrile Gloves - 12"	10	300	



Kimberly-Clark® Purple Nitrile-Xtra® Exam Gloves

Fisher Scientific Code	KCP Code	Description	Size	Total/Case	Product Features
19-149-864	55090	Kimberly-Clark® Purple Nitrile-Xtra® Exam Gloves 12" Ambi	XS	500	<ul style="list-style-type: none"> • Tested for use with 10 chemotherapy drugs and glutaraldehyde • Meet NFPA Standards 1999:2008
19-149-864A	50601	Kimberly-Clark® Purple Nitrile-Xtra® Exam Gloves 12" Ambi	S	500	
19-149-864B	50602	Kimberly-Clark® Purple Nitrile-Xtra® Exam Gloves 12" Ambi	M	500	
19-149-864C	50603	Kimberly-Clark® Purple Nitrile-Xtra® Exam Gloves 12" Ambi	L	500	
19-149-864D	50604	Kimberly-Clark® Purple Nitrile-Xtra® Exam Gloves 12" Ambi	XL	500	



Kimtech Pure® W4 Wipers – Sterile, Pre-Saturated

Fisher Scientific Code	KCP Code	Delivery System	Color	Material	Size W" x L"	Total/Case	Product Features
19-150-0300	76490	Resealable Pouch	White	Polypropylene	9" x 11"	400	<ul style="list-style-type: none"> • 100% polypropylene with 70% Isopropyl/Alcohol/30% Deionized Water • No glues or binders • Knife-cut edges minimize lint

© 2014 Thermo Fisher Scientific Inc. All rights reserved.
Trademarks used are owned as indicated at www.fishersci.com/trademarks.

In the United States:

For customer service, call 1-800-766-7000
To fax an order, use 1-800-926-1166
To order online: www.fishersci.com

