

ON YOUR TEAM

High-quality, reliable cytology solutions for your lab.

Produce consistent and accurate results
using StatLab's high quality line of cytology
stains, reagents and accessories.

Don't be locked into brand-name liquid-based cytology consumables.

Meet the StatLab cytology portfolio for gyn, urine, and FNA specimens.

Why Switch?

Utilizes common collection protocols that providers are already using for liquid-based cytology.

May be used for out of the vial testing once your laboratory has validated the StatPrep™ collection system for use with these tests.

Appropriate for use with fine-needle aspiration, urine, and gynecologic specimens.

Cost savings provide relief in an environment of declining reimbursements.

Saccomanno Cytology Collection Fluid

A ready-to-use fixative used for the collection of sputum. Lyses red blood cells. Mucoproteins remain suspended. Cells from samples are easily concentrated by centrifuging or filtering.

Saccomanno Collection Fluid, 120mL (40mL fill), 96/case STLSL26

Saccomanno Collection Fluid, 1 Gallon, 1/each STLSL251


Reagent Alcohol

50% Alcohol Prefill, 120mL, 96/case ST6950120ML

CytoVial, 95% Alcohol Prefilled Vials, 120mL, 96/case SCVP9500100

StatLab Cytology Prefills and Fixatives

Convenient for the preservation and transport of cytology specimens intended for liquid-based cytology processing. Designed to work together with your existing liquid-based cytology processor and testing to produce consistent, reliable results.



StatPrep Preserve Solutions

Used for the preservation and transport of cytology specimens intended for liquid-based cytology processing.

StatPrep Preserve Gyn 20mL, 100/case CY50GYN20ML

StatPrep Preserve Non-Gyn 20mL, 100/case CY50NONGY20

StatPrep Preserve Non-Gyn 32oz, 1/each CY32NONGY32

StatPrep Cytolysis Solution

Used for the preservation and transport of cytology specimens intended for liquid-based cytology processing.

StatPrep Cytolysis Solution, 30mL Centrifuge, 100/case CY50LYSIS30

StatPrep Cytolysis Solution, 32oz, 1/each CY32LYSIS32

CytoStat-Red Fixative

An alcohol-based general purpose fixative solution. Used for diagnostic non-gynecological cytology specimens. This fixative lyses red blood cells and solubilizes proteins while preserving the diagnostically relevant cells.

CytoStat-Red Prefill, 120mL, 96/case STLCR120ML

CytoStat-Red Prefill, 1 Gallon, 1/each STLCR1

CellBlock Fixative

Designed for cells from body fluids. Ideal for fine needle aspirations, body cavity fluids, washings, brushings, and gynecologic and non-gynecologic liquid-based specimens.

CellBlock, 1000mL 1/each STLCBF32

Cytology Fixative Spray

Economical, environmentally safe, and effective cell preservation. Natural ethyl alcohol and polyethylene glycol formulation.

Cytology Fixative Spray, 4oz, 4/case STLCYCF4CS

ON YOUR TEAM

StatPrep Liquid Based Cytology Filters

Made in a clean-room production environment with ISO quality assurance. Quality testing on every filter. Compatible with standard LBC processors. 100/case.

StatPrep Filter, Gyn	CYFILTGY100
StatPrep Filter, Non-Gyn	CYFLNONG100
StatPrep Filter, Urine	CYFILTUR100



Cytology Adhesion Slides

StatPrep Adhesion Cyto Slides, Large Circle	CYSLIDELGCI
144 slides (1 gross)	
StatPrep Adhesion Cyto Slides, Small Circle	CYSLIDESMCI
144 slides (1 gross)	
HistoBond CytoSpin White Slides, 1 Circle	STLCSP1C
1440 slides (10 gross)	
HistoBond CytoSpin White Slides, 2 Circles	STLCSP2C
1440 slides (10 gross)	

**StatPrep
Small Circle
Cyto Slide**
17mm Circle



**StatPrep
Large Circle
Cyto Slide**
33mm Circle



TissueGuard™ Processing Gel

Tissue Guard encapsulates and retains specimens during histological processing while allowing chemical reagents to penetrate the specimen. Appropriate for cytology and histology applications, TissueGuard is an alternate to other agarose stabilizers.

Box of 12 Tubes	STLTG12
-----------------	---------

Traditional Cytology Stains

StatLab Cytology stains are formulated from top-grade, certified raw materials that have been tested to meet the industry's highest standards. Finished-good testing to ensure quality, longevity and lot-to-lot consistency. All stain products are ready-to-use and require no dilution.

Hema-Diff Rapid Stain Kit, 16oz Includes Fixative, Xanthene, Thiazine	STL1197016
Hema-Diff #1 Fixative, 1 pint	STL1198516
Hema-Diff #1 Fixative, 1 gallon	STL119851
Hema-Diff #2 Xanthene, 1 pint	STL1197516
Hema-Diff #2 Xanthene, 1 gallon	STL119751
Hema-Diff #3 Thiazine, 1 pint	STL1198016
Hema-Diff #3 Thiazine, 1 gallon	STL119801
EA-65 Cytology Stain, 500mL	STLSL302
EA-50 Classic Cytology Stain, 1 pint	STLSL911
EA-50 Classic Cytology Stain, 1 gallon	STLSL9116
EA-50 Modified Cytology Stain, 1 pint	STLSL9616
EA-50 Modified Cytology Stain, 1 gallon	STLSL961
OG-6 Modified Cytology Stain, 1 pint	STLSL9216
OG-6 Modified Cytology Stain, 1 gallon	STLSL921
Reserve Cyto-DX Stain, 500mL	STLSL304



Urine Cytology Kit

Convenient method for collection, preservation, and transport of urine cytology specimens. For urine cytological testing, kit contains prefilled fixative jar, sterile 120 mL urine collection jar, and a biohazard transport bag.

Urine Cytology Kit	STLUCKTP
--------------------	----------

Pap Smear Cytology Kit

Kit includes fixative, one frosted slide, one 8" cervical brush, one dual-tip wooden scraper, pump spray fix-ative pack and detailed instructions for smear collection. Structural support protects slide during transport. 500 kits/case.

No Touch® pap kit	STLPK6000
-------------------	-----------



ON YOUR TEAM

Cytology Funnels

Fully compatible with the Thermo Shandon CytoSpin® Centrifuge. Used to deposit a thin layer of cells in a clearly defined area of a microscope slide.

Single Cytology Funnel, 500/case	STLM96410FW
Double Cytology Funnel, 500/case	STLM96420FW
Single Cytology Funnel with Filter Card, Cap and Disposable Clip, 480/case	STLM96510FW



Double Cytology Funnel



Single Cytology Funnel



Cyt-O-Dot Pens

Durable, acrylic fiber nib is beveled to a precise, fine point. The ink line or 'dot' applied is the ideal size for microscopic visualization. Permanent, waterproof, and write on almost any surface including glass, plastic, ceramic and metal. 4 pens/pack.

Cyt-O-Dot Pen Pack (Green, Blue, Red, Black)	STL9001
Cyt-O-Dot Pen Pack, Green	STL9002
Cyt-O-Dot Pen Pack, Blue	STL9003
Cyt-O-Dot Pen Pack, Red	STL9004
Cyt-O-Dot Pen Pack, Black	STL9005

Centrifuge Tubes

Translucent polypropylene. Screw-on caps provide firm, leak-resistant fit.

Centrifuge Tube, 50mL, 200/case	STLLX4103
Centrifuge Tube, 50mL, 500/case	STLLX4108
Centrifuge Tube, 15mL, 500/case	STLLX4115

HistoCool Cassette Cooling Trays

Effectively aids in the cooling of tissue blocks. Trays are filled with freezer gel that stays frozen longer: up to 6-8 hours. Simply re-freeze at the end of each day. 1/each.

Red Small Cassette Cooling Tray, 8.5x8x3in	STLSLHC20
Blue Large Cassette Cooling Tray, 12.5x8.25x4in	STLSLHC30



Cervical Broom



Scored StatPrep Cervical Brush

Sterile Cervical Broom

Cervical Broom Cell Collector, 100/pack	STLCB400100
---	-------------

Scored StatPrep Cervical Brush

Scored Cervical Brush, Sterile 8", 100/pack	STLCB500100
---	-------------

Plastic Cervical Scraper 8"

Plastic Cervical Scraper 8", 500/case	STLCT1760
---------------------------------------	-----------

Disposable Transfer Pipettes

Disposable, practical and economical. 3 mL, low density, polyethylene. Non-sterile. 500/case

Disposable Transfer Pipettes, 3mL	STLCT2000
-----------------------------------	-----------

Distributed by Fisher Healthcare. Contact us today:

In the United States

Order online: fisherhealthcare.com

Call customer service: 1-800-640-0640