

ON YOUR TEAM

# High-quality, reliable products for your entire anatomic pathology lab

---

Supporting you with consistent  
availability of the high-quality  
products your lab needs.

Find it at [fisherhealthcare.com](https://fisherhealthcare.com)

# Designed to support the highest quality of patient care.



## Leak-Proof Formalin and Specialty Prefills

High quality prefills with a StatClick™ seal prevent leaks during transit. Available in 20, 40, 60, 90 and 120 mL sizes, with 10% NBF, Alcoholic formalin, Zinc formalin, and Carson Millonig. Larger sizes available for 10% NBF prefills.

### 10% Neutral Buffered Formalin, 96/case

20mL	STLNB0507
40mL	STLNB0115
60mL	STLNB0230
90mL	STLNB0345
120mL	STLNB0460

### Buffered Zinc Formalin, 96/case

20mL	STLBFZ0507
40mL	STLBFZ0115
60mL	STLBFZ0230
120mL	STLBFZ0460

### Carson Millonig, 96/case

20mL	STLCM1507
40mL	STLCM1115
60mL	STLCM1230
120mL	STLCM1460

### NBF 5% Alcohol, 96/case

20mL	STLANB0507
40mL	STLANB0115
60mL	STLANB0230

*Additional formulas and unit quantities available including bulk configurations*

## Reagents, Decals, and Softeners

StatLab's high-performance reagents help efficiently prepare and process any type of tissue with confidence, from alcohols and xylene substitutes to decalcifiers and more.

Alcohol 100%, 1 gallon	STL69001
Alcohol 95%, 1 gallon	STL70951
Alcohol 80%, 1 gallon	STL69801
Alcohol 70%, 1 gallon	STL70701
Alcohol 50%, 1 gallon	STL70501
Decal™ Decalcifier, 32oz	STL121132
Decal Stat™ Decalcifier, 32oz	STL121232
Citra-Clear™ Xylene Substitute, 1 gallon	STL64001
XS-3™ Xylene Substitute, 1 gallon	STL74001
Methanol 99%, 1 gallon	STL403801
Isopropyl 99%, 1 gallon	STL364201



## Neutralizers

Safely neutralize and dispose of formalin while avoiding the cost of aldehyde waste recycling. StatLab's eco-friendly neutralizing products are the sensible choice for your lab and the environment.

NeutraGreen Liquid Neutralizer, 5 gallon	STL152518
NeutraGreen Liquid Neutralizer, 2.5 gallon	STL152525
NeutraGreen Liquid Neutralizer, 1 gallon	STL15251G
NeutraGreen Powder 500mL, 12/case	STLNS55012



**ON YOUR TEAM**

## Paraffin

StatLab manufactures a variety of paraffins that deliver improved infiltration and embedding, leading to more efficient deparaffinization and the minimization of wrinkles at a low melting point.

ParaPro 360, 10/case	STLPP360
ParaPro SEM, 10/case	STLPPSEM
ParaPro Blue, 8/case	STLPPBLU
ParaPro XLT, 10/case	STLPPXLT
PolyShield Paraffin Repellent 120mL, 6/case	STLPS6



## Coverslipping Tape and Coverglass

Millennia Coverslipping Tape is compatible with any specimen and is made from a proprietary, ultra clear tri-acetate film that has optical characteristics comparable to that of glass. StatLab's high-quality coverglass is especially well-suited for use in automatic coverslippers.

Millennia Coverslipping Tape	STLM7500
Coverglass 24x40mm, 10 oz/case	STLSL102440
Coverglass 24x50mm, 10 oz/case	STLSL102450
Coverglass 24x60mm, 10 oz/case	STLSL102460

*Additional sizes and unit quantities available*


## Millennia Command™ Adhesion Slides

StatLab Millennia Command™ adhesion slides help prevent tissue loss and optimize your workflow. Strong adhesion meets the rigorous demands of harsh IHC processing. Millennia Command is a hybrid slide that is suitable for IHC, special stains and H&E with no background staining. Compatible with any thermal transfer slide printer.

### Frosted 90° Corners, 100/box

White	□	STLMCOMMW90
Pink	□	STLMCOMMP90
Green	□	STLMCOMMG90
Blue	□	STLMCOMMB90
Yellow	□	STLMCOMMY90
Aqua	□	STLMCOMMA90
Lilac	□	STLMCOMML90
Tan	□	STLMCOMMT90
Red	□	STLMCOMMR90

*Additional adhesion and non-adhesion and economy slides available: ColorView Adhesion & Non-Adhesion; Millennia 1000, 2000, 2.0 Adhesion.*





## Embedding Mediums

An aqueous-based frozen section embedding medium designed to support tissue blocks in cryostat sectioning. Formulated to promote rapid freezing, enhanced sectioning and consistent results

POLARSTAT PLUS 4oz, 12/case	STLPSP12
POLARSTAT, 4oz	STLCRYO4

## Mounting Mediums

Rapid drying mounting medium, Acrymount Plus™, is compatible with xylene and xylene-substitute clearing agents including d-Limonene.

AcrymountPlus, 16oz	SL80PLUS16
AcrymountPlus, 4oz	SL80PLUS4

## FreezeSpray

StatLab FreezeSpray is ozone safe and is used to quickly cool paraffin blocks and fresh tissue in frozen sections. Two versions are available, flammable and non-flammable.

StatFreeze, 9oz	STLSL202
StatFreeze Non Flammable, 9oz	STLSL202NF

## BioMesh Screen Cassettes

New fine mesh cassettes for automatic cassette printers are ideal for small biopsies and fragile specimens, eliminating the need for messy sponges, bags, and paper. Rounded BioMesh edges ensure no tissue fragments are lost in the corners.

White Bulk Individual Cassettes, 250/pack	STLHS500
White Sleeved Cassettes, 750/case	STLHS500SL
White Threaded Cassettes, 1000/case	STLHS500T

*Additional cassette colors available*



## Supa Mega Cassettes

Ideal for embedding and processing larger tissue samples. Designed with a large open pore area to allow better reagent flow and improve turnaround times. They also ensure greater adhesion, eliminating the possibility of the embedded specimen when sectioning fibrous specimens. 100 cassettes/pack.

White Cassettes	22222015
White Mothership Cassettes	22222122

*Additional cassette colors available*

ON YOUR TEAM



## Blades

All StatLab microtome blades are developed to be ultra-durable and long-lasting, for fewer changes in the lab, leading to increased efficiency in cost and time. Blades are made in Japan with a proprietary surface ionization process, so the blade edge is slowly refined, sharpened and optimized in order to maximize hardness, sharpness and durability.

TruCut Low Profile, 50/pack	STLCUT7500
TruCut High Profile, 50/pack	STLCUT7550
StatCut Low Profile, 50/pack	STLCUT7900
StatCut High Profile, 50/pack	STLCUT7950

## Stains

A wide variety of stains for every lab's preference. StatLab manufactures standard and unique stains to the industry that quickly deliver consistent, cost-effective, high quality results to support rapid diagnosis. An accessible Technical Support team is standing by to help your customers with conversion and optimization according to their pathologist's preferences.

### Standard Stains

Hematoxylin Gill 1 Stain, 1 pint	STLSL9716
Hematoxylin Gill 2 Stain, 1 pint	STLSL9416
Hematoxylin Gill 3 Stain, 1 pint	STLSL9516
Hematoxylin Harris, 1 pint	STLSL9016
Treosin Stain, 1 pint	STLSL9316
Bluing Reagent Stain, 1 gallon	STLSL102
High-Def Acid Rinse 3.5% Solution, 1 gallon	STLSL103

### Premium/Hybrid Stains

Hematoxylin Vintage Stain, 500mL	STLSL100
Eosin-Y Vintage Stain, 500mL	STLSL101

*Additional stains and quantities available*



## TissueGuard™ Processing Gel

Tissue Guard encapsulates and retains specimens during histological processing while allowing chemical reagents to penetrate the specimen. Appropriate for cytology and histology applications, TissueGuard is an alternate to other agarose stabilizers.

Box of 12 Tubes	STLTG12
-----------------	---------





## Tissue Marking Dyes

Get the outstanding performance you expect by using StatLab marking dyes to isolate specific areas of tissue samples and identify surgical margins.

Black, 2oz bottle	■	STLSL662BK2
Blue, 2oz bottle	■	STLSL662BL2
Green, 2oz bottle	■	STLSL662GR2
Red, 2oz bottle	■	STLSL662RD2
Yellow, 2oz bottle	■	STLSL662YL2
Orange, 2oz bottle	■	STLSL662OR2
Violet, 2oz bottle	■	STLSL662VT2
Colorbond Mordant Spray 8oz 6/CS		STLCBM6

## Transport & Storage

Expertly designed for safe transport, StatLab biopsy bags are an easy choice. Efficient, cost-effective storage for slides and blocks, to help manage compliance with long-term storage guidelines.

Biopsy Bag 30X45 Small, 1000/pack	22222170
Biopsy Bag 45X75 Medium, 1000/pack	22222171
Tissue Block Storage, 25/case	22222039
Filo Slide Storage, 4 Slide Wallet, 100/pack	22222052
Filo Slide Storage, 6 Slide Wallet, 10/pack	22222057
Eco Slide Storage, 1000 Slide Storage	22222038



## Vapor Monitoring

StatLab makes it easy to monitor vapor levels of formaldehyde and xylene in your laboratory.

SafeTech Formaldehyde Monitor, 6/pack	STLVMBFORM6
SafeTech Xylene Monitor, 6/pack	STLVMBXYLN6

ON YOUR TEAM

## Traditional Cytology Stains

StatLab Cytology stains are formulated from top-grade, certified raw materials that have been tested to meet the industry's highest standards. Finished-good testing to ensure quality, longevity and lot-to-lot consistency. All stain products are ready-to-use and require no dilution.

Hema-Diff #1 Fixative, 1 pint	STL1198516
Hema-Diff #2 Xanthene, 1 pint	STL1197516
Hema-Diff #3 Thiazine, 1 pint	STL1198016
EA-65 Cytology Stain, 500mL	STLSL302
EA-50 Classic Cytology Stain, 1 pint	STLSL9116
EA-50 Modified Cytology Stain, 1 pint	STLSL9616
OG-6 Modified Cytology Stain, 1 pint	STLSL9216

*Additional cytology products and quantities available*



## Urine Cytology Kit

Convenient method for collection, preservation, and transport of urine cytology specimens. For urine cytological testing, kit contains pre-filled fixative jar, sterile 120 mL urine collection jar, and a biohazard transport bag.

Urine Cytology Kit

STLUCKTP

# StatLab has a history of reliable products, preventing disruption and allowing your team to spend more time on patient care.

0

BACKORDERS OF HEMATOXYLIN OR PARAFFIN IN 2022

### An agile, responsive partner

During the global hematoxylin and paraffin shortages in 2022, not a single StatLab customer was backordered due to raw material shortages. The experienced strategic purchasing team at StatLab observed disruptions on the horizon and took action. Working with contacts at the raw material origin, they secured additional allotments to support existing customers as well as labs impacted by shortages with other suppliers.



### Support to relieve stretched labs

25% of laboratories report staffing shortages, and 13% of histotechs in the industry now will be retiring within the next 5 years<sup>1</sup>. An experienced StatLab Technical Support team helps you with product testing, conversion and product troubleshooting, so your team can focus on supporting accurate diagnoses.

1. The U.S. Anatomic Pathology Market: Forecast & Trends 2019-21, Laboratory Economics, LLC, (2019)





Distributed by Fisher Healthcare. Contact us today:

**In the United States**

Order online: [fisherhealthcare.com](https://fisherhealthcare.com)

Call customer service: 1-800-640-0640

© 2023 Thermo Fisher Scientific Inc. All rights reserved.  
Trademarks used are owned as indicated at [fishersci.com/trademarks](https://fishersci.com/trademarks).  
BN235993352-US

