



**SPEEDY
DELIVERY.**

**SUPERIOR
SEPARATIONS.**



BUCHI FlashPure EcoFlex Silica Cartridges are now available through the Fisher Scientific channel. Small common cartridges are available for same day shipment.

Economically purify samples on your automated chromatography system without sacrificing quality.

Catalog No.	Silica Weight	Pack of	Loading Capacity	Default Flow Rate	Max. Pressure
05410003	4 g	20 pc	4.0 mg - 40 mg	15 mL/min	300 psi
05410004	12 g	18 pc	12 mg - 1.2 g	30 mL/min	300 psi
05410005	25 g	12 pc	25 mg - 2.5 g	32 mL/min	300 psi
05410006	40 g	12 pc	40 mg - 4.0 g	40 mL/min	300 psi

EcoFlex cartridges are also available to ship from BUCHI same or next day when ordered through the Fisher Scientific channel.

- 50 µm silica particles for routine separations
- Large selection of cartridge types and sizes from 4 g to 5000 g
- Non-RFID cartridges are suitable on any system from BUCHI, ISCO, Biotage, Interchim or Yamazen.

Contact us today:

In the United States

Order online: fishersci.com

Fax an order: 1-800-926-1166

Call customer service: 1-800-766-7000

In Canada

Order online: fishersci.ca

Fax an order: 1-800-463-2996

Call customer service: 1-800-234-7437

