

Contec® PROSAT® Wipes with CyQuanol™

Presaturated meltblown polypropylene wipes

PROSAT® Wipes with CyQuanol™ (meltblown polypropylene) are intended for life science cleanroom manufacturers who want to effectively maintain cGMP compliance by maximizing efficiency and minimizing inconsistency in their contamination control process.

PROSAT Wipes with CyQuanol provide a consistent release of solvent to thoroughly disinfect surfaces in critical environments. They are exceptionally clean and also free from additives of any kind. Meltblown polypropylene contains very low levels of sodium and other ions.

These wipes are EPA registered with proven efficacy claims against bacteria, viruses, and fungi. With a short contact time of 1 minute, these wipes make compliance to disinfection protocols more manageable, as microbial kill claims can be achieved without reapplication.

Using presaturated wipes delivers a controlled amount of disinfectant to surfaces every time compared to traditional

spray methods. Presaturated wipes ensure repeatability and reproducibility of disinfectant protocols to improve process consistency.

PROSAT Wipes with CyQuanol are also available with Polynit Heatseal. Both substrates are available sterile for the most critical areas.

For best results, after dwell time is achieved, wipe the surface with an IPA presaturated wipe such as PROSAT® Pi.



Features	Benefits
Meltblown polypropylene fabric	<ul style="list-style-type: none"> Exceptionally clean and free from additives
1-minute disinfectant contact time	<ul style="list-style-type: none"> Maintains compliance and appropriate dwell times
Low in particles and fibers	<ul style="list-style-type: none"> Allows for the use in critical applications
Presaturated with disinfectant	<ul style="list-style-type: none"> Reduces VOC emissions and solvent use
Resealable pouch	<ul style="list-style-type: none"> Preserves cleanliness and solvent saturation levels
Available gamma irradiated to a sterility assurance level (SAL) of 10 ⁻⁶	<ul style="list-style-type: none"> Suitable for Grade A/B (ISO 5) environments

Product	Quantity	Cat. No.
PROSAT Wipes with CyQuanol, Meltblown Polypropylene, 9" x 11" (230 x 280 mm)	30/pouch; 50 pouches/case	17-005-037
PROSAT Sterile Wipes with CyQuanol, Meltblown Polypropylene, 9" x 11" (230 x 280 mm)	30/pouch; 48 pouches/case	17-005-038

Distributed by Fisher Scientific. Contact us today:

In the United States

Order online: fishersci.com

Call customer service: 1-800-766-7000

