



# Fisherbrand™ Traceable™ has your monitoring needs covered

## Simple, digital reporting

Designed to evolve how you monitor data.

### MIN/MAX Memory

Manually monitor readings with minimum/maximum memory

### USB Data Retrieval

Download data from device to USB flash drive

### WI-FI/Ethernet Connected

Real-time, remote monitoring—receive alerts on mobile device or email



## MIN/MAX Monitoring

- Minimum/Maximum Memory
- High Low Alarms

## USB Data Logging

- Log Data Continuously
- Download Memory Data to PC

## Cloud Service Monitoring

- 24/7 Remote Monitoring
- Built in VFC Reporting

	Traceable® Min/Max Memory	Traceable® Excursion-Trac™ USB Data Logging	TraceableLIVE® Cloud Service (Wi-Fi or Ethernet)
<b>DATA STORAGE</b>	Minimum/Maximum (until cleared)	2 years (resettable)	Unlimited cloud storage
<b>USER INTERFACE</b>	Local display	Local display	Local display, Web, smartphone
<b>ALARMS</b>	Audible, visual	Audible, visual	Email, SMS text, Push notification, audible, visual
<b>DATA OUTPUT</b>	None (pen & paper)	USB data stick to Excel file	Real-time graphs, reports, alarm history, statistics, audit trail, cal due notice
<b>EMAIL/TEXT ALERTS</b>	✗	✗	✓
<b>NIST-TRACEABLE CALIBRATION CERT</b>	✓	✓	✓
<b>VACCINE FOR CHILDREN COMPLIANT</b>	✗	✓	✓
<b>21 CFR 11 COMPLIANT OPTIONS</b>	✗	✗ (available on Memory-Loc options -see fishersci.com)	✓
<b>APPLICATION</b>	<b>Cat. No.</b>	<b>Cat. No.</b>	<b>Cat. No.</b>
<b>FRIDGE/FREEZER</b>	15-079-709	15-081-123	Wi-Fi: 15-079-615   Ethernet: 15-079-704
<b>ULTRA-LOW FREEZER -80°C</b>	15-081-108	15-079-624	Wi-Fi: 15-079-626   Ethernet: 15-079-705
<b>LIQUID NITROGEN DEWAR</b>	✗	15-079-669	15-079-639
<b>ROOM MONITORING &amp; INCUBATORS</b>			
<b>TEMPERATURE</b>	15-079-708	15-081-123	15-079-616
<b>TEMPERATURE/HUMIDITY</b>	15-079-630	15-079-665	15-079-628
<b>TEMPERATURE/HUMIDITY/PRESSURE</b>	14-650-118	15-079-666	15-079-664
<b>CO<sub>2</sub>/INCUBATOR</b>	✗	✗	15-079-662

### In the United States

Order online: fishersci.com  
 Fax an order: 1-800-926-1166  
 Call customer service: 1-800-766-7000

### In Canada

Order online: fishersci.ca  
 Fax an order: 1-800-463-2996  
 Call customer service: 1-800-234-7437

