



## Fisherbrand™ Disposable Plastic Syringes

New sterile single-use, disposable syringes are individually wrapped in a paper/plastic blister package. They're available in a wide range of capacities, have a concentric tip, and offer either Luer-lock or Luer-slip tips.

- Designed for small-volume filtration and aspiration of liquids
- Suitable for use with syringe filters
- Three-part construction includes a polypropylene barrel and piston and a rubber gasket
- For research purposes only

Cat. No.	Capacity	Luer Tip	Quantity/Pack	Quantity/Case
14-955-450	1mL	Luer Slip	100	4,000
14-955-451	3mL	Luer Slip	100	3,000
14-955-452	5mL	Luer Slip	100	2,000
14-955-453	10mL	Luer Slip	100	1,200
14-955-454	20mL	Luer Slip	50	900
14-955-455	60mL	Luer Slip	50	400
14-955-456	1mL	Luer Lock	100	4,000
14-955-457	3mL	Luer Lock	100	3000
14-955-458	5mL	Luer Lock	100	2,000
14-955-459	10mL	Luer Lock	100	1,200
14-955-460	20mL	Luer Lock	50	900
14-955-461	60mL	Luer Lock	50	400

**In the United States:**

For customer service, call 1-800-766-7000

To fax an order, use 1-800-926-1166

To order online: [fishersci.com](http://fishersci.com)