



Fisherbrand[®]
QUALITY. RELIABILITY. VALUE.

Fisherbrand* Printed Disposable Culture Tubes

Ideal for tissue culture procedures


Fisherbrand Printed Disposable Culture Tubes

- Gamma-irradiated and tested for sterility
- Printed graduations
- White writing area
- Round-bottom tubes made of polypropylene
- Blue snap cap made of polyethylene
- Snap caps can be positioned for aerobic or anaerobic applications
- 14mL volume
- Three-year shelf life



Cat. No.	O.D. x L	Cap Supplied	Packaging	Quantity
14-956-9A	17 x 100mm	No	Tray of 125	1,000/Cs.
14-956-9B	17 x 100mm	Snap	Individually wrapped	500/Cs.
14-956-9C	17 x 100mm	Snap	Tray of 25	500/Cs.

Cat. No.	Description	Packaging	Quantity
14-956-9G	Sterile Polyethylene Snap Cap	Bag of 500	2,000/Cs.



14-956-9G

*Fisherband is a trademark of Fisher Scientific.



14-956-9B



14-956-9C



14-956-9A



For customer service, call 1-800-766-7000.
To fax an order, use 1-800-926-1166.
To order online: www.fishersci.com

©2012 Thermo Fisher Scientific Inc. All rights reserved.

12_0732

JD/MM-8/12

BN0830122