

## Which glove is best for you?

Kimtech™ Exam Gloves provide protection from known and unknown risks while delivering comfort, dexterity and the assurance to do your best science. They work hand-in-hand with other head-to-toe PPE solutions from Kimberly-Clark Professional. After use, give your gloves new life through The RightCycle™ Program, the first large-scale recycling program for non-hazardous PPE waste.



**Highest Protection**  
against chemotherapy drugs and other high-risk applications

LEVEL OF PROTECTION



### Kimtech™ Polaris™ Nitrile Exam Gloves

- Tested for use with chemotherapy drugs, opioid fentanyl and gastric acid (per ASTM D6978-05)
- 9.5" and 12" length
- 5.9 mil thickness for 9.5" / 7.48 mil thickness for 12" (at fingertip)
- Superior quality with AQL of .65\*



### Kimtech™ Purple Nitrile™ Exam Gloves

- Tested against chemotherapy drugs (per ASTM D6978-05)
- 9.5" and 12" length
- 5.9 mil thickness (at fingertip)
- AQL of 1.0\*



### Kimtech™ Vista™ Nitrile Exam Gloves

- Tested against chemotherapy drugs (per ASTM D6878-05)
- Approved for food contact
- 9.5" length
- 4.3 mil thickness (at fingertip)
- AQL of 1.5\*



### Sterling™ Nitrile Exam Gloves

- Tested against chemotherapy drugs (per ASTM D6978-05)
- Approved for food contact
- 9.5" and 12" length
- 3.5 mil thickness (at fingertip)
- AQL of 1.0\*



### Kimtech™ Element™ Nitrile Exam Gloves

- Exceptional value and comfort
- Approved for food contact
- 9.3" length
- 3.2 mil thickness (at fingertip)
- AQL of 1.0\*












**Basic Protection**  
against light chemical splash and other low-risk applications

\*exceeds ASTM 2.5 standard for critical defects or pinholes

# All Kimtech™ Exam Gloves are:



- FDA cleared as examination gloves
- Powder-free
- Beaded cuff for donning
- Ambidextrous
- Nitrile gloves that are not made with Natural Rubber Latex
- Recyclable through The RightCycle™ Program

Product Family	Size	Manufacturer Catalog Number	Fisher Scientific Catalog Number	Quantity per Box (10 boxes/case)	Price per Box
 <b>Kimtech™ Polaris™ Nitrile Exam Gloves</b> 9.5" Length 	XS	62770	17-103-196A	100	
	S	62771	17-103-196B	100	
	M	62772	17-103-196C	100	
	L	62773	17-103-196D	100	
	XL	62774	17-103-196E	90	
 <b>Kimtech™ Polaris™ Nitrile-Xtra™ Exam Gloves</b> 12" Length 	XS	62760	17-103-191A	50	
	S	62761	17-103-191B	50	
	M	62762	17-103-191C	50	
	L	62763	17-103-191D	50	
	XL	62764	17-103-191E	50	
 <b>Kimtech™ Purple Nitrile™ Exam Gloves</b> 9.5" Length	XS	55080	19-149-863	100	
	S	55081	19-149-863A	100	
	M	55082	19-149-863B	100	
	L	55083	19-149-863C	100	
 <b>Kimtech™ Purple Nitrile-Xtra™ Exam Gloves</b> 12" Length	XL	55084	19-149-863D	90	
	XS	55090	19-149-864	50	
	S	50601	19-149-864A	50	
	M	50602	19-149-864B	50	
 <b>Kimtech™ Vista™ Nitrile Exam Gloves</b> 9.5" Length	L	50603	19-149-864C	50	
	XL	50604	19-149-864D	50	
	XS	62825	17-000-577	200	
	S	62826	17-000-578	200	
 <b>Sterling™ Nitrile Exam Gloves</b> 9.5" Length	M	62827	17-000-579	200	
	L	62828	17-000-580	200	
	XL	62829	17-000-581	180	
	XS	50705	19-048-131	200	
 <b>Sterling™ Nitrile-Xtra™ Exam Gloves</b> 12" Length	S	50706	19-048-132	200	
	M	50707	19-048-133	200	
	L	50708	19-048-134	200	
	XL	50709	19-048-135	170	
 <b>Kimtech™ Element™ Nitrile Exam Gloves</b> 9.0" Length	XS	53137	19-050-590	100	
	S	53138	19-050-591	100	
	M	53139	19-050-592	100	
	L	53140	19-050-593	100	
	XL	53141	19-050-595	100	
 <b>Kimtech™ Element™ Nitrile Exam Gloves</b> 9.0" Length	XS	62870	17-000-582	250	
	S	62871	17-000-583	250	
	M	62872	17-000-584	250	
	L	62873	17-000-585	250	
	XL	62874	17-000-586	230	

## Our Guarantee

Your total satisfaction means everything to us. If, for any reason, our products do not meet your expectations, Kimberly-Clark will reimburse you\* for your initial purchase, via FREE product, for up to \$1,000. For more information on Kimberly-Clark Professional, visit us online at [www.KCProfessional.com](http://www.KCProfessional.com).

\*Guarantee extended to consuming end-user accounts only.

Kimberly-Clark warrants that its products (1) comply with K-C's standard specifications as of the delivery date to K-C's authorized distributors/direct purchasers (2) comply with all K-C labeling representations; and (3) are manufactured in compliance with all applicable federal, state, and local laws in effect at the time and place of manufacture of the products.

THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. K-C is not liable for any kind of special, incidental, or consequential damages. K-C's liability for breach of contract, tort or other cause of action shall not exceed the product purchase price. Purchasers and users are deemed to have accepted the above warranty and limitation of liability, and cannot change the terms by verbal agreement or by any writing not signed by K-C. To the extent required by applicable law, K-C does not limit its liability for death/injury resulting from K-C's negligence.

\*TM Trademarks of Kimberly-Clark Worldwide, Inc. or its affiliates. Marques de commerce de Kimberly-Clark Worldwide, Inc. ou de ses sociétés affiliées. © KCWW.

Distributed by Fisher Scientific. Contact us today:

In the United States

Order online: [fishersci.com](http://fishersci.com)

Call customer service: 1-800-766-7000

In Canada

Order online: [fishersci.ca](http://fishersci.ca)

Call customer service: 1-800-234-7437

